



Manual

Version 1.2

Product name: DNA Fragment Size Analysis Kit

Cat #: DFSAK-100, DFSAK-200

Description:

MCLab's DNA Fragment Size Analysis Kit offers fast and reliable separation, sizing and quantification of DNA or RNA samples utilizing Capillary Electrophoresis technology, which can separate molecules with a single base difference. The kit provides high resolution sizing and quality control of your DNA/RNA samples.

Features:

- Easy and fast sample preparation: fluorescent labeled DNA oligos can be attached to your amplified NGS library within 5 minutes.
- High resolution: separate molecules with a single base difference.
- Ultra high sensitivity: one signal for each kind of DNA molecule. Peak sizes are true to molecule quantity, no normalization is needed.
- Ultra high throughput: continuously analyze 960 samples* within a run.
*(total sample numbers depends on the Genetic Analyzer model)
- High accuracy: DNA ladder added to each sample to eliminate equipment errors and increase sizing accuracy.

Recommended Storage Condition: -20 °C

Protocol:

Fluorescent Tagging:

1. Set up the fluorescent tagging reaction in a PCR microtube:

Reagent	Volume (µl)
DNA sample (PCR product)	1-2*
Fluorescent Tag Mix	4
Total volume	5-6

2. Place the PCR microtube in a thermo cycler with the lid heating to 100°C and incubate using the following program:
 - 5min @ 37°C
 - 5min @ 96°C
 - Keep @ 4°C until next step (or place the microtube on ice)

Prepare for Capillary Electrophoresis:

1. Set up the capillary electrophoresis mix in a 96 well plate:

Reagent	Volume (µl)
Fluorescent tagging rxn	1
Red Dye standard (premixed in Super-DI)	15.5
Total volume	16.5

2. Cover the plate with a 96-well cap mat, mix the reagent vigorously on a shaker, and briefly centrifuge the plate to collect all the liquid from the side of each well.
3. Place the covered 96 well plate in a thermo cycler with the lid heating to 100°C and incubate using the following program:
5 min @ 96°C (to denature the double helix and prepare for the CE)
Keep @ 20°C until next step
4. Load the prepared 96 well plate onto a genetic analyzer, and your DNA sample is ready to be validated.

Note

- * Prepared DNA sample (after PCR) normally contains about 100ng/ul DNA fragments.
- 1. Keep both supplied reagents and the fluorescent tagging reaction in the dark at all time.
- 2. Store both supplied reagents at -20°C, thaw them on ice after removing from -20°C storage.
Keep the reagents on ice after thawing.