



# Manual

Version 2.82

**Product name:** Topomize DNA LT Library Prep Kit

**Cat #:** TOPO-100A, TOPO-100B

## The Kit Includes:

End Repair Mix 1	(Yellow)
End Repair Mix 2	(Yellow)
Topomized Adapters (12)	(Blue)
Topomization Mix	(Blue)
2x MCAmp Mix	(Red)
Primer Mix	(Red)
MCMag Beads (5 mL)	
Resuspension Buffer	(Clear)
Sequencing Primers (200 $\mu$ M)	(Green)
<u>(Read 1 / i7 Read / Read 2)</u>	

## Required Materials Not Included:

Sheared DNA sample  
Ultrasonicator for DNA shearing  
80% Ethanol (freshly prepared)  
Nuclease-free water  
Pipets and tips 10 $\mu$ L, 100 $\mu$ L, 1000 $\mu$ L  
PCR microtubes  
Microcentrifuge tubes, 1.5mL  
Microtube and microcentrifuge tube racks  
Magnetic rack/stand  
Thermocycler  
Centrifuges  
Rotator for bead mixing (optional)  
Bioanalyzer (optional)  
Horizontal electrophoresis system (optional)  
Real-time PCR system (optional)  
Fluorescence spectroscopy (optional)

## Recommended Storage Condition:

-20°C

This document and its contents are proprietary to Molecular Cloning Laboratories, Inc. (MCLAB) and its affiliates, and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of MCLAB. MCLAB does not convey any license under its patent, trademark, copyright, or commonlaw rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

MCLAB DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN OR ANY USE OF SUCH PRODUCT(S) OUTSIDE THE SCOPE OF THE EXPRESS WRITTEN LICENSES OR PERMISSIONS GRANTED BY MCLAB IN CONNECTION WITH CUSTOMER'S ACQUISITION OF SUCH PRODUCT(S).

**FOR RESEARCH USE ONLY**

© 2015-2016 MCLAB. All rights reserved.


---

## Introduction:

MCLAB's Topomize DNA Library Prep Kit is an innovative, fast and high quality DNA library construction kit for next generation sequencing. Instead of using traditional ligase-based methods, our kit uses topoisomerase-based technology to attach adapters to fragmented DNA. The Topomize DNA Library Prep Kit has unparalleled efficiency and yield, with no adapter dimers or chimeras. The kit is designed for 10ng to 1µg of sheared genomic DNA input, and is compatible with the Illumina platforms. The streamlined workflow allows the DNA library to be ready in about 90 minutes.

The purpose of this protocol is to add adapter sequences onto the ends of DNA fragments to generate indexed libraries for single-read or paired-end sequencing on the Illumina platforms.

## Workflow Process:



---

## Protocol:

**Starting material:** 10 ng-1 µg fragmented DNA

Warm up MCMag Beads to room temperature at least 30 minutes prior to use.

### 1 End Repair Reaction 1

- Mix the following components in a nuclease-free PCR microtube.

End Repair Mix 1	20 µL
DNA (fragmented) 10ng-1µg	50 µL
<hr/>	
Total	70 µL

- Mix gently by pipetting, followed by a quick spin to collect all liquid from the sides of the tube.
- Place in a thermocycler with the heated lid on, and run the following program:  
25 minutes @ 50°C  
10 minutes @ 75°C  
Hold at 4°C

### 2 End Repair Reaction 2

- Add End Repair Mix 2 directly to the previous reaction mixture.

End Repair Reaction 1	70 µL
End Repair Mix 2	8 µL
<hr/>	
Total	78 µL

- Mix gently by pipetting, followed by a quick spin to collect all liquid from the sides of the tube.
- Place in a thermocycler with the heated lid on, and run the following program:  
5 minutes @ 72°C  
Hold at 4°C
- Proceed immediately to the next step.

### 3 Topomization Reaction

- Perform the following steps on ice.
- Add the following components directly to the previous reaction mixture.

End Repair Reaction Mixture	78 µL
Topomization Mix	20 µL
Topomized adapter	2 µL
<hr/>	
Total	100 µL

- Mix gently by pipetting, followed by a quick spin to collect all liquid from the sides of the tube.
- Place in a thermocycler with the heated lid on, and run the following program:  
15 minutes @ 16°C  
Hold at 4°C
- Proceed immediately to the next step.

---

## 4 Size selection and cleanup

*(For input less than 50 ng, size selection is not recommended. Use 1x beads to perform regular cleanup.)*

- Transfer the reaction to a 1.5ml microcentrifuge tube.
- Vortex MCMag Beads until well-dispersed.
- Dilute MCMag Beads with PCR grade water to 160  $\mu\text{L}$  per 100  $\mu\text{L}$  of end-repaired sample with the following dilution formulas.

Insert size	Formula		10% excess for multiple samples	Your Calculation
350 bp	Beads	# of samples X 90 $\mu\text{L}$	# of samples X 99 $\mu\text{L}$	
	H <sub>2</sub> O	# of samples X 70 $\mu\text{L}$	# of samples X 77 $\mu\text{L}$	
550 bp	Beads	# of samples X 70 $\mu\text{L}$	# of samples X 77 $\mu\text{L}$	
	H <sub>2</sub> O	# of samples X 90 $\mu\text{L}$	# of samples X 99 $\mu\text{L}$	

- Vortex the diluted MCMag Beads until well-dispersed. Add 160  $\mu\text{L}$  diluted MCMag Beads to each tube containing 100  $\mu\text{L}$  of previous reaction from step 3, and mix thoroughly by pipetting up and down.
- Incubate at room temperature for 5 minutes, then place the tube on a magnetic stand and wait until the liquid is clear (~5 minutes).
- Transfer 250  $\mu\text{L}$  of the supernatant to a new tube (**Caution: do not discard the supernatant**).
- Vortex the undiluted MCMag Beads until well-dispersed. Add 30  $\mu\text{L}$  of the undiluted MCMag Beads to each tube, and then mix thoroughly by pipetting up and down.
- Incubate at room temperature for 5 minutes, then place the tube on a magnetic stand and wait until the liquid is clear (~5 minutes). Discard the supernatant.
- With the tube remaining on the magnetic stand, wash with 200  $\mu\text{L}$  of freshly made 80% ethanol without disturbing the beads. Pipet the ethanol out and repeat this step one more time.
- Use a 20  $\mu\text{L}$  pipette to remove the residual ethanol from each tube.
- Allow the beads to air-dry at room temperature (~5 minutes). Avoid over-drying the beads as this may result in lower recovery of target DNA.
- Add 25  $\mu\text{L}$  Resuspension Buffer to the dried beads. Remove the tube from the magnetic stand, and then mix thoroughly by pipetting. Incubate for 2 minutes at room temperature.
- Place the sample onto a magnetic stand and wait until the liquid is clear (2–5 minutes).
- Carefully transfer 23  $\mu\text{L}$  of the supernatant to a PCR microtube.



**The sample can be stored at -20°C for up to 7 days.**

---

## 5 PCR Amplification

### 5.1 Reaction

Saved Supernatant	23 $\mu$ L
2x MCAmp Mix	25 $\mu$ L
Primer Mix	2 $\mu$ L
Total	50 $\mu$ L

- Mix the above reagents gently in a PCR microtube by pipetting, followed by a quick spin to collect all liquid from the sides of the tube.
- Set up the following PCR program (with a heated lid):

	Denaturation	Amplification			Extension	Hold
Temperature	95°C	98°C	60°C	72°C	72°C	10°C
Time	10 min	10 sec	30 sec	30 sec	2 min	
		X (N) cycles				

Starting Material	(N) cycles
10 ng	10-12
50 ng	8-10
100 ng	6-8
250 ng	4-6
500 ng	3-5
1 $\mu$ g	2-4

### 5.2 Cleanup

PCR product	50 $\mu$ L
MCMag Beads	50 $\mu$ L
Total	100 $\mu$ L

- Transfer the 50  $\mu$ L post-PCR sample into a 1.5mL microcentrifuge tube and add 50  $\mu$ L of MCMag Beads. Gently pipet up and down 10 times to mix thoroughly.
- Incubate at room temperature for 5 minutes. Then place the tube onto a magnetic stand for 5 minutes, or until the supernatant is clear, and discard the supernatant.
- With the tube remaining on the magnetic stand, wash with 200  $\mu$ L of freshly made 80% ethanol without disturbing the beads. Pipet the ethanol out and repeat this step one more time.
- Use a 20  $\mu$ L pipette to remove the residual ethanol from each tube.
- Allow the beads to air-dry at room temperature for about 5 minutes. Avoid over-drying the beads as this may result in lower recovery of target DNA.
- Add 27  $\mu$ L Resuspension Buffer to the dried beads. Remove the tube from the magnetic stand, and then mix thoroughly by pipetting. Incubate for 2 minutes at room temperature.
- Place the sample onto a magnetic stand and wait until the liquid is clear (2–5 minutes).
- Carefully transfer 25  $\mu$ L of the supernatant (purified DNA library) into a new tube, without disturbing the beads.

## 6 Validation and Quantification

- The purified sample can be checked on an agarose gel or Bioanalyzer. Dilute the sample at a ratio of 1:2 with water (e.g. 3  $\mu$ L of sample with 6  $\mu$ L water) for high sensitivity Bioanalyzer chips.
- The concentration of the purified library can be checked using qPCR, with the MCNext™ SYBR® Fast qPCR Library Quantification Kit.

### Kit Contents:

- The Topomize DNA LT Library Prep Kit is available in a Set A and a Set B. Each Topomize DNA LT Library Prep Kit contains enough reagents to prepare up to 24 samples. When used together, sets A and B allow for pooling up to 24 samples using the 12 different indexes in each kit.
- The Topomize DNA LT Library Prep Kit contains 2 boxes: a Set A or Set B box and a MCMag Beads box.

#### Box 1 of 2: 24 Samples Set A or Set B Box, Store at -25°C to -15°C

Quantity	Set A Reagent	Set B Reagent	Volume
1	End Repair Mix 1	End Repair Mix 1	0.55 mL
1	End Repair Mix 2	End Repair Mix 2	0.22 mL
1	Topomization Mix	Topomization Mix	0.55 mL
1	Primer Mix	Primer Mix	60 $\mu$ L
1	2X MCAmp Mix	2X MCAmp Mix	0.7 mL
1	Resuspension Buffer	Resuspension Buffer	1.5 mL
1	Topomized adapter Index 2 (CGATGT)	Topomized adapter Index 1 (ATCACG)	5 $\mu$ L
1	Topomized adapter Index 4 (TGACCA)	Topomized adapter Index 3 (TTAGGC)	5 $\mu$ L
1	Topomized adapter Index 5 (ACAGTG)	Topomized adapter Index 8 (ACTTGA)	5 $\mu$ L
1	Topomized adapter Index 6 (GCCAAT)	Topomized adapter Index 9 (GATCAG)	5 $\mu$ L
1	Topomized adapter Index 7 (CAGATC)	Topomized adapter Index 10 (TAGCTT)	5 $\mu$ L
1	Topomized adapter Index 12 (CTTGTA)	Topomized adapter Index 11 (GGCTAC)	5 $\mu$ L
1	Topomized adapter Index 13 (AGTCAA)	Topomized adapter Index 20 (GTGGCC)	5 $\mu$ L
1	Topomized adapter Index 14 (AGTTCC)	Topomized adapter Index 21 (GTTTCG)	5 $\mu$ L
1	Topomized adapter Index 15 (ATGTCA)	Topomized adapter Index 22 (CGTACG)	5 $\mu$ L
1	Topomized adapter Index 16 (CCGTCC)	Topomized adapter Index 23 (GAGTGG)	5 $\mu$ L
1	Topomized adapter Index 18 (GTCCGC)	Topomized adapter Index 25 (ACTGAT)	5 $\mu$ L
1	Topomized adapter Index 19 (GTGAAA)	Topomized adapter Index 27 (ATTCCT)	5 $\mu$ L
1	Sequencing Primer Read 1	Sequencing Primer Read 1	20 $\mu$ L
1	Sequencing Primer Index Read	Sequencing Primer Index Read	20 $\mu$ L
1	Sequencing Primer Read 2	Sequencing Primer Read 2	20 $\mu$ L

#### Box 2 of 2: MCMag Beads Box, Store at 2°C to 8°C

Quantity	Set A Reagent	Set B Reagent	Volume
1	MCMag Beads	MCMag Beads	5 mL