

User Manual

Version 1.0

Product name: dCas13a

Cat #: dCas13a-100, dCas13a-200, dCas13a-300

Description:

The dCas13a is a mutant version of the Cas13a protein. Unlike wild-type Cas13a, it can bind to target RNA without cleaving it, making it ideal for applications where RNA localization, capture, or stabilization is required. Additionally, biotinylated dCas13a includes an N-terminal biotin tag, enabling efficient use in streptavidin-based experiments for RNA-protein interaction studies, RNA localization assays, and high-throughput screening.

Protocols:

Applications and Sample Protocol Overview

A. RNA Binding Assay

1. Design RNA guide specific to your target RNA.
2. Assemble dCas13a-guide RNA complex:
 - Mix dCas13a (1 μ g) with guide RNA (final concentration: 100 nM).
 - Incubate at 37°C for 15 minutes.
3. Introduce the complex into your experimental setup, such as live cells or lysates.

B. Streptavidin-Based Pull-Down Experiment

1. Prepare RNA sample and biotinylated dCas13a.
2. Incubate RNA sample with dCas13a in binding buffer for 30 minutes at room temperature.
3. Add streptavidin-coated beads and incubate for 1 hour at 4°C with gentle rotation.
4. Wash beads 3 times with wash buffer.
5. Elute bound complexes for downstream analysis.

C. RNA Imaging

1. Conjugate dCas13a with fluorophores using a streptavidin-fluorophore complex.
2. Use in live or fixed cells for RNA localization studies.

Troubleshooting:

| Problem | Possible Cause | Solution |
|---------------------------|--|---|
| Weak RNA binding signal | Insufficient dCas13a-guide complex formation | Optimize guide RNA concentration or sequence. |
| Loss of activity | Improper storage conditions | Use freshly thawed aliquots. Avoid freeze-thaw. |
| High background in assays | Non-specific binding of biotin tag | Add blocking steps with streptavidin. |
| Low yield in pull-down | Insufficient incubation with beads | Increase incubation time or bead concentration. |

References:

1. Abudayyeh, O. O., Gootenberg, J. S., Essletzbichler, P., Han, S., Joung, J., Belanto, J. J., ... & Zhang, F. (2017). RNA targeting with CRISPR-Cas13. *Nature*, 550(7675), 280–284. <https://doi.org/10.1038/nature24049>
2. Cox, D. B. T., Gootenberg, J. S., Abudayyeh, O. O., Franklin, B., Kellner, M. J., Joung, J., & Zhang, F. (2017). RNA editing with CRISPR-Cas13. *Science*, 358(6366), 1019–1027. <https://doi.org/10.1126/science.aaq0180>
3. Wessels, H. H., Méndez-Mancilla, A., Guo, X., Legut, M., Daniloski, Z., & Sanjana, N. E. (2020). Massively parallel Cas13 screens reveal principles for guide RNA design. *Nature Biotechnology*, 38(6), 722–727. <https://doi.org/10.1038/s41587-020-0456-9>