

# User Manual

Version 1.0  
Revision Date: 08/07/2013

**Product name:** Deoxynucleotide Triphosphate (dNTPs)

**Cat #:** dNTP-2DA, dNTP-2DC, dNTP-2DG, dNTP-2DT, dNTP-5DA, dNTP-5DC, dNTP-5DG, dNTP-5DT, dNTP-10M, dNTP-25M, dNTP-2DN, dNTP-5DN

## Description:

High quality deoxynucleotide triphosphate (dNTPs) are functionally tested in long PCR to be PCR qualified, and meet or exceed the criteria for high-quality sequencing with Thermo Sequenase DNA polymerase. High purity: 98% triphosphate, (as confirmed by chromatographic analysis and <sup>31</sup>P nuclear magnetic resonance (NMR) spectroscopy) ensures high levels of incorporation. It's supplied in buffer-free, ready to use solution and is available as individual products, sets or premixes.