

User Manual

Version 2.0
Revision Date: 09/21/2015

Product name: Terrific Broth Agar Plates

Cat #: TBAP-100, TBAP-200, TBAC-101, TBAC-102

Description:

1.2% Bacto tryptone, 2.4% yeast extract, 1.6% glycerol, 1.5% agar, 72mM K_2HPO_4 , and 17mM KH_2PO_4 .

Application:

Carbenicillin is an analog of ampicillin. For use with ampicillin resistant strains and strains harboring plasmids such as pBluescript, pGEM, pUC series plasmids.

Recommended storage condition: 4°C